






GIF Quantum™

Model 963 / 965 / 966

The GIF Quantum™ series represents Gatan's 4th and most advanced generation of post-column energy filters. Redesigned from the ground up, the GIF Quantum™ combines advanced dodecapole-based electron optics with a blazing fast CCD camera system to yield an imaging filter that defines the new state-of-the-art in the capture of both highly detailed EELS and EFTEM data sets with maximum throughput.

With the GIF Quantum™, there is no compromise between EFTEM and EELS performance. Gatan's patented sensor readout architecture (US Patents 5946033 & 7157720) allows the same CCD sensor to be used interchangeably as a full-frame high-quality imaging device, a high-speed live viewing device, and an ultra fast spectroscopy device. Aberration correction up to 5th order allows the use of a 1.8x larger 9 mm entrance aperture for EFTEM and a 2x larger 5 mm entrance aperture for energy loss spectroscopy at performance specifications superior to that of the GIF Tridium.

The GIF Quantum™ comes in several models to optimally couple to the broad range of TEM configurations and applications available today. Whether you have a work-horse LaB₆ instrument for capturing elemental maps and contrast enhanced images day after day, an extreme brightness aberration corrected STEM for probing interfaces at the atomic level, or a monochromated system for observing energy transitions at the meV level, there is a GIF Quantum™ imaging filter to fit your needs and all those in between. Advanced options such as a 1 us electrostatic shutter, a 10x faster camera, a 2 keV fast spectrum offset module, and DualEELS™ hardware and software, allow each Quantum model to be tailored to the applications at hand.



Features	Benefits
Electrostatic Shutter	Acquire images and spectra with unprecedented exposure control and dynamic range
1000 spectra per second	High-speed, dose-efficient STEM EELS spectrum imaging for detail rich mapping
9 mm entrance aperture	Large field-of-view energy filtered imaging, mapping and diffraction with narrow slit widths
2.5 & 5 mm entrance apertures	Improved collection efficiency for aberration corrected STEM EELS
Dual-speed camera readout	Low-noise imaging and high-speed viewing from the same detector
60 kV to 300 kV operation	Broad range of operating modes and configurations
2000 eV EELS range	Capture a broader range of edges in a single spectrum for simplified quantification
DualEELS™ capability	Effortlessly measure precise energy shifts and apply advanced quantification routines for a new level of EELS analysis
Advanced auto tuning	Confidence your imaging filter is operating at peak performance
Dodecapole-based optics	Outstanding energy resolution and very low image distortions

Please note: Not all GIF Quantum™ models support all features and specification. Please see the individual product datasheets for details. Specifications and features are subject to change.

Top: DualEELS™. Dual energy range EELS spectrum image data set of precipitates in steel. Data courtesy Prof. G. Kothleitner, FELMI, Graz University of Technology.
Bottom: Multi-Map. Automated multiple element mapping from a semiconductor device. 3-window maps from Si-K, Al-K, O-K, Ti-L₂₃ & N-K are automatically acquired together with a thickness map to form a complete data set.

Comparison to GIF Tridiem

60 kV minimum operating voltage	20 kV lower
0.75% maximum image distortion	2x lower
2 keV fast spectrum offset	2x larger
2 keV spectroscopy field of view	2x larger
30 fps (512 x 512) preview	6x faster
5 mm aperture for routine EELS	6¼x higher collection efficiency
2.5 mm aperture for high resolution EELS	6¼x higher collection efficiency
1,000 spectra per second	30x faster
1 µs electrostatic shutter	50,000x faster

Please note: Not all GIF Quantum™ models support all features and specification. Please see the individual product datasheets for details. Specifications and features are subject to change.

Ordering information

Part Number	Description
963	GIF Quantum™SE
965	GIF Quantum™ER
966	GIF Quantum™ERS

Please contact your Gatan sales representative for complete ordering information.

GIF Quantum™ product line

Feature	GIF Quantum™SE	GIF Quantum™ER	GIF Quantum™ERS
Energy resolution (eV FWHM)	0.25	0.10	0.04
Isochromaticity (eV P-P)	2.0	2.0	2.0
Distortion (% max)	1.5	0.75	0.75
Chromatic dist. (% 50 eV)	1.0	0.50	0.50
Aperture size (mm)	2.5 / 5	2.5 / 5 / 9	2.5 / 5 / 9
Energy field-of-view (eV)	2000	2000	2000
Spectral rate (1/s)	100	1000	1000
Fast camera modes	Optional	Included	Included
1 µs electrostatic shutter	Optional	Included	Included
Advanced BF/DF detector	Optional	Included	Included

Note: Performance at 200 kV. Specifications are subject to change.

GIF Quantum, and DualEELS are trademarks of Gatan, Inc. Gatan and the Gatan logo are registered trademarks of Gatan, Inc. Other brand and product names are the trademarks or registered trademarks of their respective owners and manufacturers.

Primary applications

- Materials research
- Failure analysis
- Catalyst research



we get it!

www.gatan.com



Certified Quality Management System

**Gatan Inc. Corporate Office/
Western USA Sales**
5794 W. Las Positas Blvd.
Pleasanton, CA 94588 USA
tel +1.925.463.0200
fax +1.925.463.0204
info@gatan.com

Eastern USA Sales
780 Commonwealth Dr.
Warrendale, PA 15086 USA
tel +1.724.776.5260
fax +1.724.776.3360
info@gatan.com

Gatan UK
25 Nuffield Way
Abingdon Oxon
OX14 1RL United Kingdom
tel. +44.0.1235.540160
fax. +44.0.1235.540169
ukinfo@gatan.com

Gatan GmbH
Ingolstädterstr. 12
D-80807 München Germany
tel +49 89 358084-0
fax +49 89 358084-77
mfelsmann@gatan.com

Gatan France
3bis, Chemin du Haut Breuil
78113 Grandchamp France
tel +33.1.34944407
GSM +33.6.80135139
fax +33.1.34871668
dmonville@gatan.com

Nippon Gatan
3F Sakurai bldg. 2-8-19
Fukagawa, Koto-ku Tokyo
135-0033 Japan
tel +81.3.5639.2772
fax: +81.3.5639.2763
rabara@gatan.com

Gatan Singapore
10 Eunios Rd. 8 #12-06
Singapore Post Centre
Singapore 408600
tel +1.65.6408.6230
fax +1.65.6293.3307
wchuang@gatan.com