

# Supported Configurations for 735.00 3D Visualization Software

---

The 3D Visualization Software product is supported by Gatan in the following configurations:

1) Dell 490 PC with graphics card  
EVGA 640-P2-N821-AR GeForce 8800GTS 640MB 320-bit GDDR3 PCI Express x16 HDCP Video Card  
Driver version: Forceware 162.18

Drivers can be found at: [http://www.nvidia.com/object/winxp-2k\\_archive.html](http://www.nvidia.com/object/winxp-2k_archive.html)