

Supported Configurations for 735.00 3D Visualization Software

Last updated: 30-Oct-09

The 3D Visualization Software product is supported by Gatan in the following configurations:

- 1) Dell 490 PC, Windows XP, 32-bit, with graphics card
EVGA 640-P2-N821-AR GeForce 8800GTS 640MB 320-bit GDDR3 PCI Express x16 HDCP Video Card
Driver version: Forceware 162.18

- 2) Dell 7400 PC, Windows XP, 32-bit, with graphics card
EVGA 640-P2-N821-AR GeForce 8800GTS 640MB 320-bit GDDR3 PCI Express x16 HDCP Video Card
Driver version: Forceware 162.18

- 3) Dell 7400 PC, Windows XP, 32-bit, with graphics card
EVGA 896-P3-1265-AR GeForce GTX 260 Core 216 896MB 448-bit GDDR3 PCI Express 2.0 x16 HDCP Ready SLI Supported Video Card
Driver version: GeForce 182.50

- 4) Dell 7500 PC, Windows XP, 32-bit, with graphics card
EVGA 896-P3-1265-AR GeForce GTX 260 Core 216 896MB 448-bit GDDR3 PCI Express 2.0 x16 HDCP Ready SLI Supported Video Card
Driver version: GeForce 182.50

- 5) Dell 7500 PC, Windows 7, 64-bit, with graphics card
EVGA 012-P3-1470-AR GeForce GTX 470 (Fermi) 1280MB 320-bit GDDR5 PCI Express 2.0 x16 HDCP Ready SLI Support Video Card
Driver version: GeForce 258.96

An archive of drivers can be found at: <http://www.nvidia.com>.